

ERGO TOUCH GLIDE

ErgoTouch GLIDE Product Information

ErgoTouch GLIDE's flat-back loaded cartridge technology allows a surgical team to smoothly manoeuvre an IOL into the cartridge with ease. When the IOL is safely manipulated into place, the cartridge's lock design ensures sleek and secure loading for safe lens delivery.

- Perfect for mono manual and bi-manual injection procedures
- Flat-back loading cartridge for easy lens loading and safe lens manipulation for the end user
- Safe lens delivery with a cartridge locking system

Product Details:

Shelf life: 3 years

Storage: From 5°C to 45°C

IOL Type Compatibility: Hydrophilic and Hydrophobic lenses

IOL Lens Compatibility: One Piece IOL

OVD Compatibility: Recommended viscoelastic Sodium ha 1.0 to 1.6 % and other viscoelastic concentration (up to 3%). Also compatible with HPMC.



Packaging Options:



ErgoTouch GLIDE Blister

- 1x Dispenser box contains 10x blister packs.
- 1x Blister pack contains 1x injector and 1x cartridge
- Dispenser box: 8.9cm width, 22.2cm height and 20.4cm depth
- Blister packs: 8.6cm width, 2.3cm height and 20.2cm depth
- Blister packs and dispenser boxes provide extra safety for both the injector and the cartridge, due to the rigid and strong packaging.

ADD THE UNIVERSAL FAMILY TO YOUR PORTFOLIO

Email: info@ergotouch.org

Item Code	Product Name	External Diameter	Incision size	
			Dock-In	Into the Bag
ErgoTouch GLIDE under BLISTER & Cartridge GMS				
ET2-B-22	ErgoTouch GLIDE under Blister with cartridges GMS 2.2	1.85	2.2	2.6



Wound-assisted technique
External Diameter: 1.6 - 2.4



Into-wound-assisted technique
External Diameter: 1.8 - 2.8



Into-bag-assisted technique
External Diameter: 2.0 - 3.0

*The incision size depends on the ophthalmic surgeon technique.